Medical Devices

AirFeel

IT > MANUALE ISTRUZIONI D'USO

EN > USER MANUAL

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NL > VERTALING VAN DE ORIGINELE

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Medical device for respiratory gymnastics

These operating instructions are provided for device model P0920EM-1.

INTENDED USE

Medical device for the prevention of increased accumulation of secretions and for therapy, in case of acute and chronic diseases of the lower airways. Prevention and/or treatment practices must be prescribed by professional medical personnel.

INDICATIONS FOR USE

Diseases of the lower airways associated with increased accumulation of mucus and secretions. It allows you to perform PEP (Positive Expiratory Pressure) therapy and O-PEP (Oscillatory Positive Expiratory Pressure) therapy.

CONTRAINDICATIONS

The device should not be used by persons with severe pneumothorax or those coughing up large quantities of blood (massive hemoptysis).

The following side effects may occur: Dizziness

If you notice any of these side effects, stop the treatment immediately. Contact your doctor or therapist to discuss how to proceed with the treatment.

Inform your doctor or therapist if you notice any side effects not listed in these instructions for use.

INTENDED USERS

The devices are intended to be used by professional medical personnel.

TARGET GROUP OF PATIENTS

The device may be used:

- on children aged 4 to 10.
- on children over 10 and on adults.

The frequency and duration of use are determined by professional medical personnel based on the patient's individual needs. Before use, the device requires a careful reading of the user manual. It is up to the medical personnel to assess the patient's condition and capabilities, to decide whether to subject him or her to the treatment. Please refer to professional medical personnel for the evaluation of the use of the device on particular types of patients, such as pregnant women, lactating women, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT

This device can be used in health care facilities, such as hospitals, clinics, etc., or even at home.

WARNINGS

- Only use the device as intended. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use.
- This product must only be used as described in these instructions for use.
- The time required to switch from storage to operating conditions is about 2 hours.
- Keep this manual safe for further reference.
- When the device is used in combination with another medical device, the instructions for use of said second device must also be complied with.
- The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely ruled out.
- The device does not administer drugs and/or substances.
- At the end of every use, once the device has been cleaned as per the paragraph HYGIENIC PREPARATION, store it in a dry and dust-free place.
- Risk of suffocation
- Some components of the device are so small that they can be swallowed by children; therefore keep the device out of the reach of children.
- Risk of ineffectiveness of therapy:
 - Use the device in the correct position (see instructions for use paragraph).
 - Use original Flaem spare parts, we decline all responsibility in the event non-original spare parts or accessories are used.
 - O-PEP treatment will be ineffective if the patient inhales/exhales too rapidly or has to force himself excessively to breathe.
- Risk of infection:
 - Before and after each use, follow the hygienic preparation operations. Make sure the device is not stored near other accessories or devices for different therapies (e.g. infusions).

TROUBLESHOOTING CASES AND RELATED SOLUTION

Problem	Cause	Solution
In case of PEP therapy: The respiratory effort is too heavy or too light	Incorrect setting of the hole on the ferrule	Check the correct position of the calibrated hole. If correct, consult your doctor / physiotherapist
In case of OPEP therapy: The sphere does not vibrate	The system is not assembled correctly	Check the correct assembly of the parts
	The ring nut is not positioned correctly	Position the ring nut in position "0"
	The cap is closed	Rotate the cap so that the indicator is aimed at the open padlock symbol

If, after having checked the conditions described above, the device still does not work properly, we recommend that you contact your trusted dealer or an authorized FLAEM service center closest to you. Find the list of all the Assistance Centers on the page <http://www.flaemnuova.it/it/info/assistance>

DISPOSAL

All components of the product can be disposed of with domestic waste unless this is prohibited by the disposal regulations in force in the respective member countries.

Packaging

Product box

Pouch packing bag

Product packaging bag

COMMUNICATION OF SERIOUS EVENTS

Please report any serious incidents that occur in connection with this device to the manufacturer and the competent authority of the Member State in which you reside.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious worsening of a person's state of health.

SYMBOLS PRESENT ON THE DEVICE OR ON THE PACKAGE

Medical device

Manufacturer

Date of production

CE marking ref. regulation
2017/745 EU and subsequent
updates

Free from phthalates and
bisphenol

Temperature limits

Atmospheric pressure limits

Consult the instructions for use

eIFU indicator

Model number

Lot code

Unique identifier of the device

Humidity limits

Attention

Sphere lock

Sphere unlock

TECHNICAL FEATURES

Model: P0920EM-1

Weight 0,110 Kg approx

MEDICAL DEVICES USABLE IN CONJUNCTION WITH THE DEVICE

The device can be used in combination with the following medical device:
Flaem nebulizer model RF7-1 or model RF7-2 (not supplied) for inhalation therapy and PEP / O-PEP therapy simultaneously.
Combination therapy must be prescribed by medical personnel.

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature	Between +10°C and +40°C
Relative humidity of the air	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature	Between -25°C and +70°C
Relative humidity of the air	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

DURATION

Model: P0920EM-1 The expected average life is 5 years

EQUIPMENT COMPONENTS AND INFORMATION ON MATERIALS

Refer to the picture on page 1 of this instruction manual.

The device equipment includes:	Material information
1 - Main body	Polypropylene
2 - Mouthpiece cap	
3 - Mouthpiece	
4 - Valve cap	
5 - Valve	
6 - PEP module	
7 - Ring nut	
8 - Ball	
9 - PEP knob	
10 - Handle	
11 - Nose clip	

INSTRUCTIONS FOR USE

Before and after each application, the components of the product must be carefully cleaned following the instructions in the paragraph HYGIENIC PREPARATION.

PEP therapy

In PEP therapy, the expiratory resistance during expiration varies according to the different sizes of the holes in the ring nut (7) to which numbers ranging from 1 to 5 are associated, the number 1 (small hole) corresponds to a higher expiratory resistance, instead the lowest corresponds to number 5 (large hole). The increase in expiratory resistance has the effect of stabilizing the airways and stimulating the dissolution of the mucus.


The size of the hole must be determined individually for each patient by professional medical personnel.

USING THE DEVICE IN PEP MODE

Preparation for therapy

Turn the PEP cap (9) to the symbol  (fig.1)

Make sure that the valve (5) is well assembled on the valve cap (4) and adheres perfectly to the internal surface of the cap (fig. 2).

Regulation of expiratory resistance

Lift the ring nut (7) and rotate it until the desired hole is aligned with the indicator  PEP module (6) (fig. 3).

Push the ring nut (7) downwards making sure that it fits correctly in the position key on the PEP module (6).

Turn the PEP Cap (9) clockwise until the position dot matches the symbol  (fig.4).

Carrying out the treatment

Please follow the steps below to administer the treatment:

- Before starting the treatment, ensure that all parts are tightly interconnected.
- Make sure that the exhalation resistance is adjusted correctly (see: Adjusting the exhalation resistance).
- If the exhalation resistance seems too high or too low during the therapeutic session, interrupt the session and then reset the device.
- Close the patient's nose with the clip supplied (11).

- The patient must hold the mouthpiece (3) between his or her teeth, hermetically sealing the mouthpiece with his or her lips.
- The patient must inhale slowly and deeply and exhale through the mouthpiece (3).
- The exhaled air must pass through the hole of the adjustment ring.
- The patient may cough to eliminate any secretions that have moved during the treatment.
- For reasons of hygiene, the patient must avoid coughing in the AirFeel system.

O-PEP therapy

The patient repeatedly exhales into the device against a moving ball, causing a vibration that is transmitted to the lungs. This vibration tirs the mucus in the lower airways. During use, the airways are kept open and mucus clearance is enhanced.

If Airfeel is held horizontally, the frequency will fluctuate during use and equal to approximately 15 Hz. This frequency can be varied by tilting the device up (higher frequency) or down (lower frequency). This allows to select the ideal frequency for each individual patient.

USE OF THE DEVICE IN O-PEP MODE

Preparation for therapy

Turn the PEP cap (9) to the symbol  (fig.1)


Make sure that the valve (5) is well assembled on the valve cap (4) and adheres perfectly to the internal surface of the cap (fig. 2).

Regulation of expiratory resistance

Lift the ring nut (7) and rotate it until the "0" symbol is aligned with the indicator  present on the PEP module (6). (fig.5).

Push the ring nut (7) downwards making sure that it fits correctly in the position key on the PEP module (6).

Make sure that the PEP Cap (9) is in position , if not, turn it anticlockwise until it is reached.

Execution of the treatment

Please follow the steps below to administer the treatment:

- The patient must sit in an upright position and relax, grab the handpiece (10) of AirFeel and hold it vertically to the floor.
- The patient must hold the mouthpiece (3) between his or her teeth and hermetically seal it with his or her lips.
- NOTE: Make sure that the holes of the upper part are not obstructed or covered by the hand.
- The patient must inhale as slowly and deeply as possible through the nose and exhale

slowly and fully through the mouth of the device (as slowly and for as long as possible). While he or she exhales, he or she must tighten the muscles of the cheeks and keep them tight.

- The patient must repeat the inhalation and exhalation cycle
- End the treatment session if the effort becomes excessive.

Notes for optimal treatment:

Place a free hand against the chest while the patient is inhaling and exhaling: when the device is used correctly, an unmistakable regular vibration can be felt. The intensity of the vibration varies according to the device's alignment with the floor: align the device in the position in which the vibration is felt more strongly (slightly angle the device upwards or downwards).


Therapy in combination with a medical device model RF7-1 or RF7-2

It is feasible both in the case of PEP therapy and O-PEP therapy (for all the instructions and restrictions on use of the device model RF7-1 or RF7-2, refer to the nebulizer's manual provided).

USE OF THE DEVICE IN COMBINED MODE WITH MEDICAL DEVICE MODEL RF7-1 or RF7-2

Having defined the therapy to be applied PEP or O-PEP and carried out the operations described in the respective paragraphs "Preparation for the therapy" and "Adjustment of the expiratory resistance" and "Execution of the treatment" connect the output of the nebulizer model RF7-1 or RF7-2 in the valve cap (4) as shown in the figure.

In this mode, the mouthpiece to be used is the one supplied with the AirFeel device and not the one with the exhalation valve supplied with the Flaem nebuliser model RF7-1 or RF7-2. Do not use any masks supplied with the nebuliser model RF7-1 or RF7-2.

HYGIENIC PREPARATION

Take apart the AirFeel system in its individual parts.

Sanitization

Sanitize the components by choosing one of the methods provided in the table and described below.

method A: Sanitize the accessories under hot drinking water (about 40°C) with a delicate dishwashing detergent (non-abrasive).

method B: Sanitize the accessories in the dishwasher with a hot cycle (70°C).

method C: Sanitize the accessories by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly with hot drinking water (about 40°C).

Disinfection

After having sanitized the components, disinfect them by choosing one of the methods provided in the table and described below. Each method can be used for a limited number

of times (see data in the table).

method A: Get a disinfectant of electrolytic chloroxidizing type (active ingredient: sodium hypochlorite), specific for disinfection, available in all pharmacies.

Execution:

- Fill a container of suitable size to contain all the individual components to be disinfected with a solution based on drinking water and disinfectant, respecting the proportions indicated on the disinfectant package.
- Completely immerse each single component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Retrieve the disinfected components and rinse them thoroughly with warm drinking water.
- Dispose of the solution according to the directions provided by the disinfectant manufacturer.

method B: Disinfect the components by boiling in water for 10 minutes; use demineralized or distilled water to avoid limescale deposits.

If you also want to carry out sterilisation, skip to the Sterilisation section.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Sterilisation

Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

At the end of each use, store the device in a dry and dust-free place.

Table of planned methods / patient accessories									
Accessory patient									
Method	1	2	3	4-5*	6	6-7*	7	8	9
Sanitisation									
method A	✓	✓	✓	✓	✓	\	✓	✓	✓
method B	✓	✓	✓	✓	✓	\	✓	✓	✓
method C	✓	✓	✓	✓	✓	\	✓	✓	✓
Disinfection									
method A	300	300	300	300	300	\	300	300	300
method B	300	300	300	300	300	\	300	300	300
Sterilization	300	300	300	300	\	300	\	300	300
✓: expected \: not expected 300: ✓ MAX 300 TIMES * Assembled									

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© 2025 FLAEM NUOVA* All rights reserved
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.